

## Protocol for RoosterVial™-hBM-1M-XF EXPANSION

Expansion of RoosterBio XF hBM-1M-XF to yield 10 million cells	
RoosterVial-hBM-1M-XF Lot No:	Date/Time:

### 1.0 MATERIALS

ITEM	VENDOR	PART No*	EXP DATE
RoosterVial-hBM-1M-XF, (1 million (M) cells)	RoosterBio	MSC-031	
CTS TrypLE Select Enzyme	Life Technologies	A1285901	
DPBS (without Ca <sup>++</sup> , Mg <sup>++</sup> )			
DPBS (without Ca <sup>++</sup> , Mg <sup>++</sup> ) + hPL or SPENT MEDIA for quench			
hMSC High Performance Media Kit XF	RoosterBio	KT-016	
2 x T225 CELLBIND flasks OR 6 x T75 CELLBIND flasks	Corning		

\* Vendors and part numbers are included for critical items.

### 2.0 MEDIA PREPARATION & CELL EXPANSION

- 2.1 Bring hMSC High Performance Media Kit XF to room temperature. Add 1 vial hMSC Media Booster XFM (SU-016) to 500 mL hMSC High Performance Basal Media (SU-005).
- 2.2 Obtain RoosterVial-hBM-1M-XF from liquid nitrogen dewar and immediately thaw in 37°C water bath. Monitor the process and remove from water bath once a small bit of ice is remaining (2-3 min).
- 2.3 Spray vial well with 70% isopropyl alcohol before transferring into biosafety cabinet.
- 2.4 Aseptically transfer cells into a 50 mL centrifuge tube.
- 2.5 Slowly (dropwise) add 4 mL of culture media to the cells.
- 2.6 Centrifuge at 200 x *g* for 10 min.
- 2.7 Carefully remove the supernatant without disturbing the cell pellet.
- 2.8 Resuspend the cells in 5 mL of culture media. When cells are resuspended bring volume up to 30 mL with culture media.
- 2.9 Mix well and seed cells equally into two T225 or into six T75 vessels, and add media to bring volume up to final volume:

Type of culture vessel (X)	Total volume of cell suspension per vessel	Cells/ cm <sup>2</sup>	Final Volume/ flask
T75 x 6 □		2,222	15ml
T225 x 2 □		2,222	45ml

- 2.10 Transfer vessel(s) into 37°C incubator and ensure that the surfaces are covered with media.
- 2.11 Microscopically observe culture everyday from day 3 onwards to determine percentage confluency.
- 2.12 If culture is less than 50% confluent on day 3, perform a media change. Completely remove the spent media from the vessel, and replace with same volume of fresh culture media. Transfer vessel back into incubator.
- 2.13 When culture is >80% confluent, prepare to harvest the following day.

Day	Confluency (%)
3	
4	
5	
6	
7	

### 3.0 CELL HARVESTING

- 3.1 For harvesting, transfer vessel into biosafety cabinet and remove spent media.  
\*Collect ~10 mL spent media in sterile container if using to quench harvest enzyme.
- 3.2 Remove media and add 10 mL of TrypLE to T225 or 3 mL TrypLE to each T75 flask, and incubate in 37°C incubator.
- 3.3 Check culture every 5 min until cells are detached from surface. Gently tap to dislodge remaining cells from surface.

Total time required for cells detachment

- 3.4 Add equal volume of quench or spent media to stop the TrypLE activity.
- 3.5 Transfer the cell suspension into a 50 mL centrifuge tube.

Total volume of cell suspension (=A)

- 3.6 Centrifuge at 200 x g for 10 min.
- 3.7 Aspirate the supernatant and resuspend cells with 4-5 mL of fresh media.
- 3.8 Measure the total volume of cell suspension.
- 3.9 Mix well and transfer 0.1 mL of cells into microcentrifuge tubes for cell counts.
- 3.10 Dilute cells to 0.5 mL with /DPBS to get counts in the range of 0.1-1 x 10<sup>6</sup> cells/mL.
- 3.11 Mix well and cells are ready for counts with cell counting device.

Raw data		Adjusted data	
Dilution Factor (=B)	Viable Cell Concentration (=C)	Cell concentration (D)=B*C	Total cells at harvest (E)=D*A

- 3.12 Place harvested cell suspension in 4°C refrigerator for use within 1 hour of harvest.