

# BOOST CONFIDENCE IN EV COLLECTION

## ROOSTERCOLLECT™-EV | COLLECTION MEDIA EXTRACELLULAR VESICLE (EV) PRODUCTION

LOW-PARTICULATE MEDIUM  
STREAMLINED. SCALABLE.



ROOSTERBIO.COM

## RoosterBio® Radically Simplifying Use of MSCs



A low particle medium, RoosterCollect-EV supports seamless transition from hMSC expansion to EV collection.

Catalyze a complete extracellular vesicle (EV) hyperefficient manufacturing workflow by starting with RoosterBio hMSC bioprocess systems including high-volume xeno-free hMSCs and paired media. This engineered system paired with RoosterCollect™-EV increases EV yields with minimized processing times allowing you to get to your product. Productivity translates readily from 2D (batch) to 3D (fed-batch) applications to scale with your development goals.

**We provide process recommendations to boost your EV yields.**

AS PART OF A COMPLETE SYSTEM FOR MSC EV COLLECTION, ROOSTERCOLLECT-EV IS ENGINEERED TO STREAMLINE INTEGRATION INTO YOUR EV WORKFLOW.

### Transition from hMSC expansion to EV collection




Fig 1(a)

### RoosterCollect™-EV reduced particle levels ensure confidence




Fig 1(b)

Fig 1(a). Switch from cell expansion phase, using RoosterBio's high-volume hMSCs and paired RoosterNourish-MSC-XF expansion medium, directly to collection phase using RoosterCollect-EV.

Fig 1(b). RoosterCollect-EV is a media engineered for extracellular vesicle collection with a particle level significantly lower than commercially available media - typically containing background particle noise.

# Low Particle Collection Medium for Use in Dynamic EV Manufacturing Processes

## APPLICATIONS FOR ROOSTERCOLLECT™-EV

Extracellular Vesicles | Secreted Proteins | Conditioned Media

### ROOSTERCOLLECT-EV SUPPORTS SCALABLE EV COLLECTION FROM hMSCs

#### EV collection supported over 48 hours




Fig 2(a)

#### Particle collection in EV size ranges




Fig 2(b)

Fig 2(a). RoosterCollect-EV supports EV collection from MSCs from 2 to 48 hours. Particle level in the conditioned media increases over the collection phase.

Fig 2(b). Collected particles have diameters in the size range of 50-350nm.

#### EV collection scales




Fig 2(c)

RoosterCollect-EV particle collection scales from 2D vessel culture to 3D bioreactor systems which increases EV production as higher cell quantities are achieved.

#### PRODUCT FEATURES

- > Low Particle Content
- > Xeno-free. Protein Free. Chemically Defined.
- > Streamlines EV Collection
- > Collect in 2D and 3D Applications
- > Component of EV Production

#### BENEFITS

- > Reduced Processing Times
- > Ensure End Product Purity

**SIMPLIFY EV PRODUCTION & ACCELERATE PRODUCT DEVELOPMENT**



**RoosterBio**<sup>®</sup>  
Radically Simplifying Use of MSCs

ROOSTERBIO.COM | 301.200.5366 | INFO@ROOSTERBIO.COM