



Donor-Specific Angiogenic Effects of MSCs in 3D Microvasculature Tissue Constructs

Hannah A. Strobel, Ph.D.
Staff Scientist
Advanced Solutions Life Sciences

Elie Zakhem, Ph.D.
Process Development Manager
RoosterBio

Introduction

Vascularization is an essential element for tissue homeostasis and function. Oxygen, for example, can diffuse approximately 100 - 200 μm into a tissue from an oxygen source, depending on the tissue. The vasculature is critical in bringing the oxygen source (i.e., blood) in close approximation to the cells of a tissue. Not only does a vessel network transport blood throughout thick tissues, but vascular cells and associated perivascular niche cells play a critical role in overall tissue function. Yet, tissue vascularization remains a key challenge throughout the field of tissue biofabrication and tissue modeling. An important aspect of tissue vascularization is the interplay between the mesenchymal stem/stromal cells (MSCs) of the tissue stroma and the vasculature. It's well known that MSCs condition the tissue environment with a variety of growth factors and cytokines, including angiogenic factors. Additionally, the vascular and perivascular cells of microvessels influence tissue cell behavior important in tissue function.

Increasingly, Angiomics[®] isolated human adipose-derived microvessel fragments (haMVs) are being used to vascularize 3D tissues and tissue models, including adipose, liver, and bone [1-3]. Angiomics[®] haMVs, isolated from discarded lipoaspirates, retain their native structure, microvessel phenotype, and cellular complexity [4, 5]. When embedded in a 3D matrix, haMVs undergo sprouting angiogenesis with growing neovessels arising from parent microvessel fragments forming an interconnected neovascular network [1, 5]. When provided a perfusion source, isolated microvessel-derived neovasculatures rapidly inosculate with that source and mature into a persistent microcirculation [6].

RoosterBio's human mesenchymal stem/stromal cells (hMSCs) are widely used in a variety of applications including tissue modeling and advanced therapy manufacturing [7, 8]. In addition to phenotypic plasticity, important in establishing tissue function, MSCs support tissue vascularity [9, 10], in part by promoting angiogenesis [11, 12]. RoosterBio hMSCs secrete a spectrum of pro-vascular growth factors, cytokines, and paracrine factors.

Here, we examined the potential of hMSC lots from different donors to promote angiogenesis in haMV-derived microvasculatures. We focused on specifically

Benefits of Combining Angiomics[®] haMVs with RoosterBio hMSCs to Vascularize Tissues


- Accelerated angiogenesis and neovascular growth
- Applicable to a variety of tissue types and systems
- Readily add cellular complexity to tissues
- Higher fidelity to the native, in vivo tissue environment
- Compatible with BioAssembly[®] automation and high-content imaging analyses

the paracrine activity of hMSCs from three different donors, reflecting differences in growth factor secretion profiles, on sprouting angiogenesis from the haMVs. 3D constructs of haMVs were co-cultured with hMSCs suspended above the constructs using cell culture inserts. This allowed the cells and vessels to share the medium compartment, but do not allow for physical interaction. The three hMSC lots promoted angiogenesis, the extent of which depended on the lot, highlighting the benefits of using both technologies in tissue modeling and regenerative applications.

METHODS

Xeno-Free RoosterVial[™] bone marrow derived hMSCs from RoosterBio [MSC-031] were cultured using RoosterNourish[™]-MSC-XF expansion media [part no. KT-016] according to the manufacturer's instructions. Three different lots reflecting 3 different donors were used, designated 257, 277, and 280. hMSCs were suspended in 3 mg/ml collagen type I at a density of 1 million cells per ml and dispensed into Transwell[®] inserts (200 μl per insert) with three inserts per group. While gelling, ASLS Angiomics[®] haMVs [ASLS-000123] were thawed, suspended at a density of 100k haMVs/ml in 3 mg/ml cold neutral collagen. In wells of a 24 well plate, 300 μl of cold, neutral acellular collagen was dispensed per well. Before this collagen could gel, 100 μl of the haMV/collagen suspension was dispensed into the middle of each well. This creates a central region of haMVs surrounded by a ring of acellular collagen. The inserts containing gelled hMSC/collagen constructs were suspended above wells containing the haMV constructs and

culture medium (RPMI + B27) was added to each well, such that both the hMSC insert and haMVs are covered (**Figure 1**). A control group contained haMVs only without hMSC inserts. After 6 days of culture, haMV constructs were fixed, stained with fluorescent UEA lectin, and imaged using an confocal scanner. Vessel length densities (total vessel length per tissue area) in constructs, a measure of angiogenesis, were determined using BioSegment™ (ASLS).



▲ **Fig. 1** (above) Experimental configuration for paracrine factor testing of hMSCs. hMSCs, embedded in collagen, are suspended above a construct of haMVs in collagen with both compartments sharing culture medium.


RESULTS AND DISCUSSION

In the present study, we explored the potential of hMSCs to stimulate angiogenesis from pre-existing microvessel segments in 3D in vitro constructs. We have previously observed the pro-angiogenic effects of hMSCs on angiogenesis in a tissue model [1]. Here, we expanded this study to focus specifically on paracrine-mediated angiogenesis by hMSCs from multiple donors. The cell culture insert format of this experiment physically separated the hMSCs and haMV constructs, allowing the hMSCs to condition the culture medium in the absence of direct cell-cell contact (**Figure 1**). Importantly, no exogenous growth factors were added. Thus, any effects on angiogenesis observed reflects the cell activity in the system.

After 6 days, haMV constructs with hMSC inserts exhibited many growing angiogenic sprouts originating from the parent vessels of the haMV isolate as compared to the haMV alone cultures (**Figure 2**).

Differences in vessel length density between control haMV constructs and those co-cultured with the three hMSC lots were measured from confocal images of lectin-stained cultures (**Figures 3-4**). Vessel length density, calculated by summing the length of each individual vessel in the image and dividing by the total area per image, is a common measure of angiogenic growth. Lot 257 hMSCs promoted the most angiogenesis (largest increase in vessel length density), followed by 277 and 280. While lot 257 hMSCs were most strongly angiogenic, all lots promoted angiogenesis.

Because there is no direct hMSC-microvessel contact, the pro-angiogenic effect of the hMSCs is due to paracrine factors and other signaling molecules released by hMSCs



▲ **Fig. 2** (above) Angiogenesis in an haMV construct. Parent microvessels were isolated from discarded adipose tissue and embedded in collagen matrix. As compared to the absence of hMSCs (A), angiogenesis, neovessel sprouts arising from parent microvessels, is enhanced in the presence of hMSCs (B).


Table 1. Growth factors secreted by RoosterBio hMSCs (pg/10⁵ cells/day)

Lot #	bFGF	HGF	IL-8	TIMP-1	TIMP-2	VEGF
257	14.5	882.3	151.3	98592.3	10694.5	687.4
277	27.9	508.4	239.8	107253.8	9009.1	592
280	21.5	671	226.3	185800.8	18188.2	725.3

into the culture medium. Examining the growth factor secretion profiles of the three hMSC lots used in this study indicates the cells are producing a variety of growth factors, many of which are known to stimulate angiogenesis (**Table 1**). Thus, it is not surprising that angiogenesis is stimulated in the hMSC presence.


Interestingly, Lot 257, which stimulated angiogenesis the most, did not produce the highest levels of factors, except for HGF. HGF is also a known pro-angiogenic growth factor [13], and may be contributing to the higher angiogenic potential of this lot. Of course, there are a variety of factors produced by hMSCs that were not measured, including small molecule metabolites and cytokines [14], that could be driving the angiogenesis from haMVs. Additionally, as fragments of native microvasculature, the haMVs are comprised of multiple cell types in addition to endothelial cells, including pericytes and macrophages. These other cell types may be responding to the hMSC-derived factors by releasing their own, angiogenic factors. An important aspect of this assay is that the cellular complexity reflects more the in vivo tissue environment, and all assay readouts reflect the interplay between pro- and anti-angiogenesis generating a “net angiogenic” outcome.

The Angiomics® haMVs have been used to vascularize a variety of tissues and tissue models, including adipose,



▲ **Fig. 3** (above) Vascular length density of haMV cultures without and with hMSC lots (257, 277, or 280). *P<0.05, one way ANOVA with Tukey post-hoc test, n=3.

bone, liver, endocrine, and cardiac [15-19], and we have previously demonstrated their dynamics for in vitro studies or implantation purposes [20, 21]. In all cases, endogenous or included parenchyma and stromal cells contribute to microvessel dynamics in the absence of exogenous factors. Here, we show that hMSCs similarly accelerate angiogenesis and neovessel growth from the haMVs without the need for growth factor supplementation. Thus, combining RoosterBio’s hMSCs with haMVs has broad applicability related to tissue construction and vascularization for in vitro 3D tissue modeling and potentially clinical applications as well. In addition, differences in hMSC angiogenic potential enable fine tuning/optimizing tissue vascularization strategies. Furthermore, the study highlights the utility of haMV cultures as an angiogenesis assay in which neovessel arise from parent microvessels, as it happens in the native tissue.



▲ **Fig. 4** (above) Representative confocal maximum projections of microvessel cultures conditioned with different lots of MSCs. During analysis, only intact parent and neovessels were counted towards vessel length density (arrows). Red = lectin stain (endothelial cells).

Authored by

Hannah A. Strobel, Ph.D.

Staff Scientist

Advanced Solutions Life Sciences

Elie Zakhem, Ph.D.

Process Development Manager

RoosterBio

References

1. Strobel, H.A., T. Gerton, and J.B. Hoying, Vascularized adipocyte organoid model using isolated human microvessel fragments. *Biofabrication*, 2021.
2. Moss, S.M., et al., A Biofabrication Strategy for a Custom-Shaped, Non-Synthetic Bone Graft Precursor with a Prevascularized Tissue Shell. *Front Bioeng Biotechnol*, 2022. 10: p. 838415.
3. Moss, S.M., et al., Point-of-use, automated fabrication of a 3D human liver model supplemented with human adipose microvessels. *SLAS Discov*, 2022. 27(6): p. 358-368.
4. Nunes, S.S., et al., Implanted microvessels progress through distinct neovascularization phenotypes. *Microvasc Res*, 2010. 79(1): p. 10-20.
5. Hoying, J.B., C.A. Boswell, and S.K. Williams, Angiogenic Potential of Microvessel Fragments Established in Three-Dimensional Collagen Gels. *In Vitro Cell. Dev. Biol.-Animal*, 1996. 32: p. 409-419.
6. Shepherd, B.R., et al., Rapid perfusion and network remodeling in a microvascular construct after implantation. *Arterioscler Thromb Vasc Biol*, 2004. 24(5): p. 898-904.
7. Timsina, H., et al., A comparative evaluation of layer-by-layer assembly techniques for surface modification of microcarriers used in human mesenchymal stromal cell manufacturing. *Biotechnology Journal*, 2022. 17(8): p. 2100605.
8. Visconti, R.J., K. Kolaja, and J.A. Cottrell, A functional three-dimensional microphysiological human model of myeloma bone disease. *J Bone Miner Res*, 2021. 36(10): p. 1914-1930.
9. Takebe, T., et al., Engineering of human hepatic tissue with functional vascular networks. *Organogenesis*, 2014. 10(2): p. 260-7.
10. Yamamoto, K., et al., Construction of Continuous Capillary Networks Stabilized by Pericyte-like Perivascular Cells. *Tissue Eng Part A*, 2019. 25(5-6): p. 499-510.
11. Lee, J.H., Y.S. Han, and S.H. Lee, Long-Duration Three-Dimensional Spheroid Culture Promotes Angiogenic Activities of Adipose-Derived Mesenchymal Stem Cells. *Biomol Ther (Seoul)*, 2016. 24(3): p. 260-7.
12. Bhang, S.H., et al., Transplantation of cord blood mesenchymal stem cells as spheroids enhances vascularization. *Tissue Eng Part A*, 2012. 18(19-20): p. 2138-47.
13. Morishita, R., et al., Therapeutic Angiogenesis using Hepatocyte Growth Factor (HGF). *Current Gene Therapy*, 2004. 4(2): p. 199-206.
14. Aghazadeh, Y., et al., Cell-based therapies for vascular regeneration: Past, present and future. *Pharmacol Ther*, 2021: p. 107976.
15. Moss, S.M., et al., Point-of-use, automated fabrication of a 3D human liver model supplemented with human adipose microvessels. *SLAS Discov*, 2022.
16. Strobel, H.A., T. Gerton, and J.B. Hoying, Vascularized adipocyte organoid model using isolated human microvessel fragments. *Biofabrication*, 2021. 13(3).
17. Moss, S.M., et al., A Biofabrication Strategy for a Custom-Shaped, Non-Synthetic Bone Graft Precursor with a Prevascularized Tissue Shell. *Frontiers in Bioengineering and Biotechnology*, 2022. 10.
18. Aghazadeh, Y., et al., Microvessels support engraftment and functionality of human islets and hESC-derived pancreatic progenitors in diabetes models. *Cell Stem Cell*, 2021.
19. Sun, X., et al., Transplanted microvessels improve pluripotent stem cell-derived cardiomyocyte engraftment and cardiac function after infarction in rats. *Sci Transl Med*, 2020. 12(562).
20. Nunes, S.S., et al., Angiogenic potential of microvessel fragments is independent of the tissue of origin and can be influenced by the cellular composition of the implants. *Microcirculation*, 2010. 17(7): p. 557-67.
21. Strobel, H.A., et al., Stromal Cells Promote Neovascular Invasion Across Tissue Interfaces. *Front Physiol*, 2020. 11: p. 1026.



RoosterBio[®]

1.301.200.5366

RoosterBio.com

info@RoosterBio.com



5295 Westview Drive, Suite 275

Frederick, MD 21703

Advanced
SOLUTIONS

1.877.438.2741

www.advancedsolutions.com



1901 Nelson Miller Pkwy

Louisville, KY 40223