

## Recommended Expansion Protocol for Fed-Batch Culture Regimen

### Protocol Summary

To expand one vial of xeno-free, human bone marrow or umbilical cord-derived Mesenchymal Stem/Stromal Cells (RoosterVial-hBM/hAD/hUC-1M-XF) using 2 Ambr<sup>®</sup> 250 Modular, multi-parallel bioreactors, you will need the following reagents, materials, and equipment:


(The protocol may be scaled proportionally to accommodate up to 8 Ambr<sup>®</sup> 250 bioreactor vessels.)

### Materials & Equipment

Item	Quantity	Vendor	Part Number*
RoosterVial-hBM/hAD/hUC-1M-XF	1 Vial	RoosterBio	MSC-031/C46001AD/C43001UC
RoosterNourish™-MSC-XF	2 Bottles	RoosterBio	K82016
RoosterReplenish™-MSC-XF	1 Bottle	RoosterBio	SU-023
T225 CellBIND flasks	2 Flasks	Corning	3293
Low Concentration Synthemax™ II Microcarriers	7.0 g	Corning	3781
Ambr <sup>®</sup> Vessel, Single Impeller, Unbaffled	2 unit	Sartorius	001-2A33
TrypLE™ Select Enzyme	1 Bottle	Life Technologies	12563029
DPBS (without Ca <sup>++</sup> , Mg <sup>++</sup> )	1 Bottle	Life Technologies	14190144
10 mL Costar Wide Tip Serological Pipettes	1 Pack	Corning	4492
100 µm Cell Strainer	5 Units	BD Falcon	352360
150 mL Sterile Bottle	1 Bottle	Corning	431175
500 mL Centrifuge Bottle	1 Bottle	Corning	431123
Biosafety Cabinet			
Centrifuge			
Incubator			
Water Bath (or ThawSTAR <sup>®</sup> )			

### Process Overview

#### Cell Expansion Summary



*\*RoosterBio strongly recommends the use of Corning CellBIND surfaces for expansion of the Xeno-Free product line.*

*\*\*Please refer to protocol for full process instructions.*

CONFIDENTIAL. FOR SOLE USE OF RECIPIENT. ALL PRODUCTS ARE FOR RESEARCH USE ONLY.

## Recommended Protocol

### 1. Expansion Options

Vessel	Surface Area (cm <sup>2</sup> )	Number of Vessels Needed	Total Surface Area (cm <sup>2</sup> )	Seeding Density (cells/cm <sup>2</sup> )	Approx. Yield at Harvest	Days of Culture
T225	225	2	450	2,200	≥ 12M	3 to 6

### 2. Media Preparation

- 2.1. Bring RoosterNourish-MSX components to room temperature, protected from light, for up to four hours.
- 2.2. Prepare 1 bottle of medium by aseptically adding 1 bottle of RoosterBooster™-MSX (Part No. SU-016) to 1 bottle of RoosterBasal™2.0-CC (Part No. M22520).
- 2.3. Mix well by capping and gently mixing the bottle.

### 3. Cell Thawing & Seeding

- 3.1. Aseptically transfer 10 mL of prepared medium into a 15 mL centrifuge tube
- 3.2. Thaw RoosterVial in an automated thawing device (e.g., ThawStar), or manually in a 37°C water bath. When thawing in a water bath, monitor the vial closely and remove from water bath once only a small bit of ice is remaining (2-3 min).
- 3.3. Aseptically transfer vial into a Biosafety Cabinet (BSC).
- 3.4. Transfer vial contents into the 15 mL centrifuge tube containing prepared medium and mix cell suspension well.
- 3.5. Centrifuge at 350 x g for 6 min.
- 3.6. Aspirate the supernatant and resuspend cells in 10 mL of RoosterNourish-MSX medium.
- 3.7. Mix well and seed cells equally into two T225 vessels, and add medium to bring volume up to final volume according to table below:

Type of Culture Vessel	Total Volume of Cell Suspension per Vessel	Final Volume per Vessel
T225	5 mL	45 mL

- 3.8. Transfer vessels into an incubator (37°C, 5% CO<sub>2</sub>) and ensure surfaces are covered evenly with media.

### 4. Cell Expansion

- 4.1. Microscopically monitor cell confluency starting on day 3 of culture.
- 4.2. When culture is >80% confluent, cells are ready to harvest.

Day	3	4	5	6
Cell Confluency				

**Note:** For best expansion and functional performance, it is recommended to passage the cultures before reaching 90% confluence. If the cultures reach over confluence, this may result in increased aggregation, decreased cell viability, growth inhibition and loss of differentiation potential.

### 5. Ambr® 250 Modular Setup & Media Conditioning

CONFIDENTIAL. FOR SOLE USE OF RECIPIENT. ALL PRODUCTS ARE FOR RESEARCH USE ONLY.

**Note:** Steps 5.1-5.6 should be completed one day before Section 6. **Flask Harvest & Bioreactor Inoculation**

- 5.1. Write a recipe in the Ambr® definition software. Considerations for the recipe are based on the following recommended parameters for hMSCs (RoosterVial™-hBM-XF or RoosterVial™-hUC-XF) in the Ambr® 250 bioreactor system:
  - Agitation direction: Down
  - Agitation speed: To facilitate cell attachment, 6 – 10 cycles of static (0 rpm for 45 min) and dynamic (100 rpm for 2 min) after bioreactor inoculation. The agitation speed should be maintained at 100 rpm for the proceeding 18 hours, and then increased by 25 rpm every 12 hours.
  - CO2 Flow: 5% of headspace gas flow (Initiated after DO spot calibration)
  - DO(%) Saturation: 100% (Initiated after DO spot calibration)
  - Headspace gas flow (Air / mix): 20 mL / min
  - pH: pH control is only monitored for this process
  - Temperature: 37°C
- 5.2. Open the pre-written recipe in the Ambr® runtime software. The user interface will prompt the following:
  - Check process definition
  - Confirm connected gases
  - Enter spot calibration data
  - Clean reservoir slots
- 5.3. In a BSC, unpackage bioreactor vessels and aseptically add 150 mL of RoosterBasal™2.0-CC medium into each of the 2 bioreactor vessels.
- 5.4. Transfer vessels out of the BSC and place on the Ambr® 250 Modular units. Continue the recipe to start the run and initiate media conditioning at the set parameters.
- 5.5. Perform a DO spot calibration to 100% and confirm initiation of DO controller at 100% saturation and headspace CO2 flow at 5% of headspace gas flow.
- 5.6. After 4 hours of media conditioning, perform pH calibration:
  - Take 1 mL medium sample through the sampling port.
  - Take an offline pH reading.
  - Input offset pH data into the Ambr® software.
  - Initiate pH control
- 5.7. In a BSC, weigh 3.5 g of sterile microcarriers into a sterile bottle that allows for complete collection of liquid/microcarrier suspension for each bioreactor vessel.
- 5.8. Add 90 mL of complete medium (RoosterBasal™2.0-CC + RoosterBooster-MSC-XF) to the microcarriers in each bottle and swirl to wet the microcarriers and equilibrate them for culture as per manufacturer's recommendation.
  - 5.8.1. If no manufacturer recommendations are provided, store the bottle containing microcarriers and complete medium overnight at 4°C.

## 6. Flask Harvest & Bioreactor Inoculation

*\*The following instruction reflects the steps required for each vessel and can be performed in parallel for all 8 vessels.*

- 6.1. Warm bottles containing microcarriers to room temperature.
- 6.2. Add 3ml of RoosterBooster-MSC-XF XF to each bioreactor unit containing 150 mL of RoosterBasal™2.0-CC.
  - 6.2.1. Borrow the vessels from the modular units and transfer to a BSC. Open the septum cap and aseptically add 3 ml of RoosterBooster-MSC-XF to each vessel.
  - 6.2.2. Return the vessels to the modular unit.
- 6.3. Transfer the T225 flask into a BSC and remove spent media.
- 6.4. Add 10 mL TrypLE per T225 flask.
- 6.5. Distribute TrypLE evenly to cover all the cells and place flasks in 37°C incubator. Check culture every 5 minutes until cells are detached from surface (typically 10-15 minutes).

CONFIDENTIAL. FOR SOLE USE OF RECIPIENT. ALL PRODUCTS ARE FOR RESEARCH USE ONLY.

6.6. Gently tap to dislodge remaining cells from surface.

**Total Time Required for Cell Detachment**

- 6.7. Add equivalent volume of RoosterNourish™-MSC-XF to each vessel to stop the TrypLE activity.
- 6.7.1. If the freshly harvested cell solution contains significant cell aggregates or extracellular matrix clumps due to overconfluence, filter solution through a 100 µm cell strainer to remove.
- 6.8. Transfer the cell suspension into a sterile 50 mL centrifuge tube.
- 6.9. Centrifuge at 350 x g for 6 min.
- 6.10. Aspirate the supernatant.
- 6.11. Resuspend cells in ~10 mL medium. Measure the total volume of cell suspension:

**Total Volume of Cell Suspension (=A)**

- 6.12. Transfer 0.5 mL of cell suspension into microcentrifuge tubes for cell counts.
- 6.13. Count cells with a cell counting device, performing a dilution if required to get within its acceptable range:

Raw Data		Adjusted Data	
Dilution Factor (=B)	Viable Cell Concentration (=C)	Cell Concentration (D)=B*C	Total Cells at Harvest (E)=D*A

- 6.14. For each Ambr250 bioreactor, proceed with the following steps, working quickly:
- 6.14.1. Add  $5.8 \times 10^6$  cells to each bottle that now contains the microcarriers and expansion media.
- 6.14.2. Mix contents well.
- 6.14.3. Add bottle contents to the AMBR250 vessel.
- 6.14.4. Engage vessel in the modular unit.
- 6.15. Proceed with recipe to facilitate cell attachment to microcarriers using the intermittent seeding method: 8 cycles of static culture for 45 minutes, with 2 minutes of agitation at 100 RPM in between cycles to redistribute the cells and microcarriers.
- 6.15.1. Turn off the DO controller during the static cycles to ensure no inaccurate control is initiated when the sensor is covered by the settled microcarriers.
- 6.16. Following the cell attachment phase, continue recipe with agitation speed at 100 RPM for 24 hours following inoculation, followed by a 25 RPM speed increase every 12 hours.

## 7. Bioreactor Feeding: Addition of RoosterReplenish™-MSC-XF

- 7.1. On day 3 of culture, add 5 mL of RoosterReplenish-MSC-XF (Part No. SU-023), 2% volumetric addition, to the bioreactor as a feed.
- 7.1.1. Borrow the vessels from the modular units. In a BSC, open the septum cap and aseptically add 5 ml of RoosterReplenish to each vessel.
- 7.1.2. Return the vessel to the modular unit.

## 8. Bioreactor Sampling

CONFIDENTIAL. FOR SOLE USE OF RECIPIENT. ALL PRODUCTS ARE FOR RESEARCH USE ONLY.

- 8.1. Borrow the vessels from the modular units and transfer into a biosafety cabinet.
- 8.2. Open the bioreactor septum cap and using a wide tip serological pipette aspirate 3 mL of uniform cell/microcarrier suspension into a 15 mL tube. (Additionally, 1 mL of cell/microcarrier suspension can be sampled to monitor cell attachment, growth, and aggregation via microscopy). Ensure cell / microcarrier suspension are homogeneous when sample is drawn by pipetting the cell-microcarrier suspension up and down first and sampling when the cell-microcarriers are suspended.
- 8.3. Return the bioreactor vessels to the modular unit and continue recipe.
- 8.4. For the 3 mL cell count sample, allow cells and microcarriers to settle into bottom of tube by gravity, typically 5-10 minutes.
- 8.5. Once the sample of cells and microcarriers settle, carefully remove as much spent medium/ supernatant from the tube, without disturbing the cell/microcarrier suspension.
- 8.6. Add ~2.7 mL (equivalent to the volume of media removed from the sample; the volume of cells / microcarrier is typically 0.3 mL) of TrypLE solution to conical tube with cells/microcarriers and gently mix, then incubate for 15 minutes at 37°C for cell dissociation from microcarriers. If large cell clumps remain, mix gently with pipette, and return to 37°C for an additional 10 minutes.
  - 8.6.1. For Day 1 and Day 2 Sampling, it may be necessary to concentrate the sample by reducing the volume of TrypLE so that the cell count falls within the optimal range of the cell counting device.
- 8.7. When cells are dissociated from microcarriers, separate cells from the microcarriers using a 100 µm cell strainer. Sample the cell suspension to perform a cell count to determine cell density in bioreactor.

Day	Raw Data		Adjusted Data	
	Dilution Factor (=B)	Viable Cell Concentration (=C)	Cell Concentration (D)=B*C	Total Cells at Harvest (E)=D*500 mL
3				
4				
5				
6				

**Note:** Cells are typically ready to harvest on day 5-6 of culture. Cultures should be harvested during exponential growth phase (i.e., not when the cells have reach a growth plateau).

## 9. Bioreactor Harvest

- 9.1. Transfer cell culture vessels and other necessary materials to biosafety cabinet.
- 9.2. Allow the cells/microcarriers to settle to the bottom of the bioreactor vessels.
- 9.3. Remove the bioreactor cap and aspirate as much spent media as possible from the culture without removing the cells/microcarriers.
- 9.4. Add 200 mL DPBS into the bioreactor and swirl gently to wash the cells/microcarriers.
- 9.5. Allow the cells/microcarriers to settle to the bottom of the bioreactor vessel.
- 9.6. Aspirate as much DPBS as possible from the culture without removing the cells/microcarriers.
- 9.7. Add 200mL of TrypLE into the bioreactor (>200 mL).
- 9.8. Return the vessels to the modular units and continue agitation at 400 RPM for 30 minutes at 37°C.
- 9.9. Transfer the bioreactor into the biological safety cabinet and measure the total volume of cells by transferring the cells/microcarrier solution into a sterile 500 mL bottle:

CONFIDENTIAL. FOR SOLE USE OF RECIPIENT. ALL PRODUCTS ARE FOR RESEARCH USE ONLY.

**Total Volume of Cell Suspension (=A)**

- 9.10. Pipette to break up remaining cell/microcarrier clusters.
- 9.11. Filter solution through a 100 µm cell strainer.
- 9.12. Quench TrypLE with equivalent volume of fresh media.
- 9.13. Mix cell solution and obtain final cell count for the bioreactor harvest.

Raw Data		Adjusted Data	
Dilution Factor (=B)	Viable Cell Concentration (=C)	Cell Concentration (D)=B*C	Total Cells at Harvest (E)=D*2A

- 9.14. Process with downstream processing as applicable.

*Caution to Users: RoosterBio products contain human sourced materials and should be treated as potentially infectious. Employ universal safety precautions and wear protective clothing and eyewear while handling. Practice appropriate disposal techniques per CDC guidelines for biohazardous material.*

*Provision of Seller Product subject to Seller Standard Terms and Conditions. Any technical advice furnished, or recommendation made concerning any use or application of any Seller Product is believed to be reliable, but Seller makes no warranty, either express or implied, as to its accuracy or completeness or of the results to be obtained. Purchaser assumes full responsibility for quality control, testing and determination of suitability of Seller Product for its intended application or use.*

CONFIDENTIAL. FOR SOLE USE OF RECIPIENT. ALL PRODUCTS ARE FOR RESEARCH USE ONLY.